



Reference: PTXBC063384

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	immunoglobulin heavy constant delta
GENE ID	3495
ORIGIN SPECIES	Human
ACCESSIONS	BC063384
VECTOR	pDONR223
NCBI SYMBOL	IGHD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr