

C1orf201-chromosome 1 open reading frame 201 Gene

Reference: PTXBC063891

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0490 protein C1orf201; C1orf201; MAPO2; O(6)-methylguanine-induced apoptosis 2; O6-methylguanine-induced apoptosis 2; sperm-tail PG-rich repeat-containing protein 1; sperm tail PG-rich repeat containing 1		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 201
GENE ID	90529
ORIGIN SPECIES	Human
ACCESSIONS	BC063891
VECTOR	pDONR223
NCBI SYMBOL	C1orf201
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttccctcaatgtgcgccgattggacaccatcatttctaataccctgcagcgaaatgcataactatccatcggaatttatttccaagagagactttagtaattcgtgtccagcatgttccaggttccaagctttatgaaagctcgaagtttgaaactcctgcaccaactattacaatgcctctgtctcttgcagcagagagaacaacgtctgtactcagcggtttatgtcaaaaacccaaagaggatcttgcgttttgcgtgataaaggacctccccagggcattatgatatacaacgaatccctgtgaagcagtcgcgaatacattaatgtctgttttaatacaaaaacacacgtggatataaactgacgtcaacaggcccggtgacgtgatttacaacccagtgatgcacaaagtccaataaaagactctttcccgaaaaacccatcctgaacttctgtcagccttgcctctgctcgaagccactttccaggtctgtcagatgagatcgtggactacttaggcccccgcaagcatttcatctctagtcacatcgtgtccaataccagccggtggacagcgccgcctcagccagcgctggtgcccagctacgtacaagccagagcttccaggaagcagtccttctctacaacaggacagaagaatggatcccggttctctag