

MRPL14-mitochondrial ribosomal protein L14 Gene

Reference: PTXBC065005

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L14mt; L32mt; MRP-L14; MRP-L32; MRPL32; RMPL32; RPML32; 39S ribosomal protein L14, mitochondrial; 39S ribosomal protein L32, mitochondrial; mitochondrial ribosomal protein L14		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L14
GENE ID	64928
ORIGIN SPECIES	Human
ACCESSIONS	BC065005
VECTOR	pDONR223
NCBI SYMBOL	MRPL14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctttcttactgggctctggggcccttcacctgtgtaagcagagtctgagccatcactgtttcagcaccactgggagctctgagtgcgattcagaagatgacgcgggtacagtggtggacaacagtgcctggggaacagccatccatcggtctctcgtgcatccatgtctat
aagaagaatggagtgggcaaggtggcgaccagatacttgccatcaagggaagagaagaaaggcgtcattgtggggcactgcatgctgccccgaatgacccccagattcgactccaacaacgtggtctcattgaggacaacgggaacctgtggggacacgaattaag
acacccatccccaccagcctgcgaagcgggaaggcgagtattccaaggtgctggccattgtcagaactttgtgtga