

CNPY2-canopy 2 homolog (zebrafish) Gene

Reference: PTXBC065015

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HP10390; MSAP; TMEM4; ZSIG9; protein canopy homolog 2; MIR-interacting saposin-like protein; canopy 2 homolog; transmembrane protein 4; canopy FGF signaling regulator 2		

For research use only.

GENE DESCRIPTION	canopy 2 homolog (zebrafish)
GENE ID	10330
ORIGIN SPECIES	Human
ACCESSIONS	BC065015
VECTOR	pDONR223
NCBI SYMBOL	CNPY2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaggctgggttgctggccctgtctctggggccctgtgggaaccgctgggctggaggagccagatctccactgtggagcatgcagggtctgtggatgaactagaatggaaattgccagggtggacccaagaagaccattcagatgggatcttccggatcaatcc
agatggcagccagtcagtggtggaggtgccttatgccgctcagaggccacctcacagagctgtgaggagatgtgaccggatgaaggagtatggggaacagattgatcctccaccatcgcaagaactacgtacgtgtagtggccgggaatggagaatccagtgaactggac
ctacaaggcatccgaatcgactcagataatagcggcaccctcaagtttgcgtgtgagagcattgtggggaatacagaggatgaactcattgaattctttccgagaggctgacaatgttaagacaaactttgcagtaagcgaacagatcttgtgacatgccctgcacatcgcacatga
gctatga