



Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

SERGEF-secretion regulating guanine nucleotide exchange factor Gene

Reference: PTXBC065375

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DELGEF; Gnefr; secretion-regulating guanine nucleotide exchange factor; deafness locus associated putative guanine nucleotide exchange factor; guanine nucleotide exchange factor-related protein; secretion regulating guanine nucleotide exchange factor		

For research use only.

GENE DESCRIPTION	secretion regulating guanine nucleotide exchange factor
GENE ID	26297
ORIGIN SPECIES	Human
ACCESSIONS	BC065375
VECTOR	pDONR223
NCBI SYMBOL	SERGEF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr