

CLK4-CDC-like kinase 4 Gene

Reference: PTXBC065732

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein serine threonine kinase Clk4; dual specificity protein kinase CLK4; CDC like kinase 4		

For research use only.

GENE DESCRIPTION	CDC-like kinase 4
GENE ID	57396
ORIGIN SPECIES	Human
ACCESSIONS	BC065732
VECTOR	pDONR223
NCBI SYMBOL	CLK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgccattccaaaagaactcactgtctgattgggatagcagagaagctggggacatgaaagctatcgtggaagtcacaaagcggagaggagatctcatagtagcacacaagagaacaggcattgtaaacacatcaccagtttaagaatctgattgtaa