

ZNF449-zinc finger protein 449 Gene

Reference: PTXBC065938

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZSCAN19; zinc finger protein 449; zinc finger and SCAN domain-containing protein 19		

For research use only.

GENE DESCRIPTION	zinc finger protein 449
GENE ID	203523
ORIGIN SPECIES	Human
ACCESSIONS	BC065938
VECTOR	pDONR223
NCBI SYMBOL	ZNF449
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtggccctgggtgtgcaatccaggcatcctgaatcaaggctctgtgttcaagaatgatactgactgtgaagttttccgtcagcgttcagcagttccagtacagagaagcagctgggcctcatgaagcatttaacaaactctgggagctttgtgtcaatggctgaagccaaa
gatgcgtctaaagaaacaaatcctggagctgtagtgttgagcaattcctaactatcctgccacagagatagagacctgggtgaggagcactgccagagaatagagaagagttgtgtcactgatagaagacttacagagagaacttgatataccagagcagacacatttgcctgt
gaatacagtcactactag