

FAM168B-family with sequence similarity 168, member B Gene

Reference: PTXBC066347

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM168B; MANI; myelin-associated neurite-outgrowth inhibitor; family with sequence similarity 168 member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 168, member B
GENE ID	130074
ORIGIN SPECIES	Human
ACCESSIONS	BC066347
VECTOR	pDONR223
NCBI SYMBOL	FAM168B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatcctgtttatagtcctggatcttctgggttcctatgcaaatgccaaaggaattggtatccagctggtttcccatgggctatgcagcagcagctcctctattctcctaacaatgtatcctggagcgaatcctctccaaacaggttacactcctggcacaccttacaagtgtcctgtt
ccccaccagcgggctgtgccaccgtactcctctcaccacagctgtgcccagagtcctgtgacaccagtgggcaaggcagctgggccagatctcaggccagccgtttgtctcctagcagggttctgtgctgcccacacggagagccctagagagcctcatggattgtaac
taa