

SFRS2-splicing factor, arginine/serine-rich 2 Gene

Reference: PTXBC066958

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SFRS2; PR264; SC-35; SFRS2A; SRp30b; serine/arginine-rich splicing factor 2; SR splicing factor 2; splicing component, 35 kDa; splicing factor SC35; splicing factor, arginine/serine-rich 2; serine and arginine rich splicing factor 2		

For research use only.

GENE DESCRIPTION	splicing factor, arginine/serine-rich 2
GENE ID	6427
ORIGIN SPECIES	Human
ACCESSIONS	BC066958
VECTOR	pDONR223
NCBI SYMBOL	SFRS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctacggccgccccctccgatgtggagggtatgacctcccaaggtggacaacctgacctaccgacctgccccacacgctgaggcgctcttcgagaagtacggcgcgctggcgacgtgtacatcccggggatcgctacaccaaggagtcgccggcttcgcttc
gttcgcttcacgacaagcgcgacgtgaggacgtatggatgccatggagcgagccgagtcggagccgttcaggctcgcagccgatctcgtacagccgctcgaagtctcggtcccgacactgttctcgalctcggctgacctccaagtcagatccgcacgaaggtcca
agtccaagtcctcgtcgggtccagatctcgttcggtccaggtcccggtctcgggtccagagtcctccccaggtgtccaagaggggaatccaaatccaggtcgcgatcgaagagtccccccaagtcctcgaagaggaaggagcggtgtcctcttaagaaatgctgcggtctcgtttg
a