

SCGB2A2-secretoglobin, family 2A, member 2 Gene

Reference: PTXBC067220

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGB1; UGB2; mammaglobin-A; mammaglobin 1; secretoglobin family 2A member 2		

For research use only.

GENE DESCRIPTION	secretoglobin, family 2A, member 2
GENE ID	4250
ORIGIN SPECIES	Human
ACCESSIONS	BC067220
VECTOR	pDONR223
NCBI SYMBOL	SCGB2A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattgctgatgctcctcatgctggcgccctctccagcactgtacgcaggctctggctgccccctattggagaatgtgattccaagacaatcaatccacaagtgtctaaagactgaatacaaaagaactcttcaagagttcatagacgacaatgccactacaatgccatagatgaattg
aaggaaatgttttctaaccaaacggatgaaactctgagcaatgttgagggtttatggttaatttcattttctctataagctttttaaatccccctgaccagggacaagtgaggctcttcaattctcactgacaatgccaaagccactagtgaacaagccttttctacattggttaa