



Reference: PTXBC067358

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARF GAP GIT1; ARF GTPase-activating protein GIT1; CAT-1; CAT1; G protein-coupled receptor kinase interacting ArfGAP 1; G protein-coupled receptor kinase-interactor 1; GRK-interacting protein 1; cool-associated and tyrosine-phosphorylated protein 1; GIT Ar		

GENE DESCRIPTION	G protein-coupled receptor kinase interacting ArfGAP 1
GENE ID	28964
ORIGIN SPECIES	Human
ACCESSIONS	BC067358
VECTOR	pDONR223
NCBI SYMBOL	GIT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]