



Reference: PTXBC067861

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAPB; DRT; EK5; EPHT3; ERK; Hek5; PCBC; Tyro5; ephrin type-B receptor 2; EPH-like kinase 5; developmentally-regulated Eph-related tyrosine kinase; elk-related tyrosine kinase; eph tyrosine kinase 3; protein-tyrosine kinase HEK5; renal carcinoma antigen NY		

GENE DESCRIPTION	EPH receptor B2
GENE ID	2048
ORIGIN SPECIES	Human
ACCESSIONS	BC067861
VECTOR	pDONR223
NCBI SYMBOL	EPHB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]