

LOC645683-ribosomal protein L13a pseudogene Gene

Reference: PTXBC067891

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADMIO; ADMIO1; APRF; HIES; signal transducer and activator of transcription 3; DNA-binding protein APRF; acute-phase response factor		

For research use only.

GENE DESCRIPTION	ribosomal protein L13a pseudogene
GENE ID	645683
ORIGIN SPECIES	Human
ACCESSIONS	BC067891
VECTOR	pDONR223
NCBI SYMBOL	LOC645683
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgcccaagaccagcgaggccatgcctccttgactgcctcaaggtgtttgacggcatccaccgccctacataagaaaaagcggatgggtgttctgctgcctcaaggtcgtgctctgaagcctacaagaaagtttgcctctggggcgacaggctcaagaggttagg
tggaagtaccaggcagtgacagccaccctggaggagaagaggagaagcaagatccactactggaagaaacagctcatgaggtctacgaaacaggccgagaagaacgtgaagaaaaactga