



Reference: PTXBC068533

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD62; CD62P; GMP140; GRMP; LECAM3; PADGEM; PSEL; CD62 antigen-like family member P; GMP-140; granule membrane protein 140; granule membrane protein 140kDa; granulocyte membrane protein; leukocyte-endothelial cell adhesion molecule 3; platelet activation d		

GENE DESCRIPTION	selectin P (granule membrane protein 140kDa, antigen CD62)
GENE ID	6403
ORIGIN SPECIES	Human
ACCESSIONS	BC068533
VECTOR	pDONR223
NCBI SYMBOL	SELP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]