

CRLS1-cardiolipin synthase 1 Gene

Reference: PTXBC069010

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf155; CLS; CLS1; GCD10; dJ967N21.6; cardiolipin synthase (CMP-forming); cardiolipin synthase 1		

For research use only.

GENE DESCRIPTION	cardiolipin synthase 1
GENE ID	54675
ORIGIN SPECIES	Human
ACCESSIONS	BC069010
VECTOR	pDONR223
NCBI SYMBOL	CRLS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtagccttgccgctggcgccgctgctggggccctgcgcggcgccgcttgggtccgggaacgcggccgagtaagcgaacgcctgctgggcccctgctgcgcggccgctgctgggcccctgctgctgggctgctggcgaacgctggaggctgcgtccggcgcttggcttgcg
gctgcccgggagcggccagcgggaaccactgttcgggcgagggaagcggctccagggccagcggcgagcggcgccgctgcccgaagccccggcgccagtgggggccggcgaacccccagcctgtatgaaacccatggacaatcccgaatatgtgtcaatgacg
agaattggcttggcccgatctggtgctattgattatgaagaagatttaataatgacactaggagttttgctttagctggactaacagatttggatggatttattgctcgaactgggccaatcaaatcagcttgggaagtgccttgcacttgcctgataaaactatcagatcttata
tgtagcttgacatgcagatcttatccagcgaacactgccaagtattcaatccttgcctatgccactgtag