



Reference: PTXBC069228

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPI transamidase component PIG-S; GPI transamidase subunit; phosphatidylinositol glycan, class S; phosphatidylinositol-glycan biosynthesis class S protein; phosphatidylinositol glycan anchor biosynthesis class S		

For research use only.

GENE DESCRIPTION	phosphatidylinositol glycan anchor biosynthesis, class S
GENE ID	94005
ORIGIN SPECIES	Human
ACCESSIONS	BC069228
VECTOR	pDONR223
NCBI SYMBOL	PIGS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]