

SCGB1D1-secretoglobin, family 1D, member 1 Gene

Reference: PTXBC069289

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LIPA; LPHA; LPNA; secretoglobin family 1D member 1; lipophilin A (uteroglobin family member); lipophilin-A; prostatein-like lipophilin A		

For research use only.

GENE DESCRIPTION	secretoglobin, family 1D, member 1
GENE ID	10648
ORIGIN SPECIES	Human
ACCESSIONS	BC069289
VECTOR	pDONR223
NCBI SYMBOL	SCGB1D1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctgcgggtgtctctcgtcgtcagctggcccttctgtctaccgggcaaatgcagtggtctgccaagctcttgggtctgaaatcacaggctcttattagctggaaaacctgtgtccaagttccaacttgccaaatttaaggcacctctggaagctgttcagccaagatggaagtga
agaaatgcgtggatacgatggcctatgagaaaagagtgctaattacaaaacattgggaaaaatagcagagaaatgtgatcgtga