



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

GDNF-glial cell derived neurotrophic factor Gene

Reference: PTXBC069369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HFB1-GDNF; ATF; ATF1; ATF2; HSCR3; glial cell line-derived neurotrophic factor; astrocyte-derived trophic factor; glial cell derived neurotrophic factor		

For research use only.

GENE DESCRIPTION	glial cell derived neurotrophic factor
GENE ID	2668
ORIGIN SPECIES	Human
ACCESSIONS	BC069369
VECTOR	pDONR223
NCBI SYMBOL	GDNF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr