

## CCL24-chemokine (C-C motif) ligand 24 Gene

Reference: PTXBC069391

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | Ckb-6; MPIF-2; MPIF2; SCYA24; C-C motif chemokine 24; CK-beta-6; chemokine (C-C motif) ligand 24; eosinophil chemotactic protein 2; eotaxin-2; myeloid progenitor inhibitory factor 2; small inducible cytokine subfamily A (Cys-Cys), member 24; small-inducib |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chemokine (C-C motif) ligand 24                                                                   |
| <b>GENE ID</b>          | 6369                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC069391                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CCL24                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcaggcctgatgaccatagtaaccagcctctgttcttgggtgtgtgcccaccacatccctacggcctgtgtggtcatccctctccctgctgcatgttcttgtttccaagagaattcctgagaaccgagtggtcagctaccagctgtccagcaggagcacatgcctcaaggcagg  
agtgatcttcaccaccaagaaggccagcagttctgtggcgacccaagcaggagtggtccagaggtacatgaagaacctggacgccaagcagaagaaggcttccctagggccaggcagtggtgtcaaggccctgtccagagatactgtgcaaccaaacctgctaa