

HIST1H4L-histone cluster 1, H4I Gene

Reference: PTXBC069392

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H4/m; H4FM; H4M; histone H4; H4 histone family, member M; Histone 4 family, member M; histone 1, H4i; histone cluster 1, H4i; histone family member; histone cluster 1 H4 family member i		

For research use only.

GENE DESCRIPTION	histone cluster 1, H4I
GENE ID	8368
ORIGIN SPECIES	Human
ACCESSIONS	BC069392
VECTOR	pDONR223
NCBI SYMBOL	HIST1H4L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgggcgcggaaggcgggaagggtctgggcaaggaggcgctaagcgccaccgcaagttctgcgcgacaattcaggcctacacagcccgccatccgacgcctggcagcggtggaggcgtaagcgcatctcaggccttatatacaggagacacgcggagttct
taaaagtgttttggagaatgtaatccgcgatgcagttacctacacggagcagcgaacgcaagacagtcacagccatggacgtgtttacgcgctcaagcgccaggcgccacccgtatggcctttggcggtga