



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RP11-265F14.2-elastase 2B Gene

Reference: PTXBC069412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RP11-265F14.2; ELA2B; chymotrypsin-like elastase family member 2B; elastase 2B; pancreatic elastase IIB (ELA2B); chymotrypsin like elastase family member 2B		

For research use only.

GENE DESCRIPTION	elastase 2B
GENE ID	51032
ORIGIN SPECIES	Human
ACCESSIONS	BC069412
VECTOR	pDONR223
NCBI SYMBOL	RP11-265F14.2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]