



Reference: PTXBC069416

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GENEY; TRIFIC; tripartite motif-containing protein 46; gene Y protein; tripartite, fibronectin type-III and C-terminal SPRY motif protein; tripartite motif containing 46		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 46
GENE ID	80128
ORIGIN SPECIES	Human
ACCESSIONS	BC069416
VECTOR	pDONR223
NCBI SYMBOL	TRIM46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]