



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

TM4SF5-transmembrane 4 L six family member 5 Gene

Reference: PTXBC069519

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane 4 L6 family member 5; tetraspan transmembrane protein L6H; transmembrane 4 superfamily member 5; transmembrane 4 L six family member 5		

For research use only.

GENE DESCRIPTION	transmembrane 4 L six family member 5
GENE ID	9032
ORIGIN SPECIES	Human
ACCESSIONS	BC069519
VECTOR	pDONR223
NCBI SYMBOL	TM4SF5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatgcacggaaaatgtgccctctgtgtggggctctccctcataccctctgcctctgtctgattgtggccaaacgcctctctgtgtaactaatggggagacctcctggaccaacaccaacatctacgtgtgcaagtctggctcatggggggcttattggcggggccctaagtgtactgt
gtccaggggatgtcagccgctgtgggcagggggcgaagggtctgtgtgggtctgtgggtgtgtggaaacccgttcaggagatgtctgcctcgtgtcttctctggcgttcgggggtctgtgtggccattactgcctctcgggtctgtgagctggggctccgaatggaccacagatgcttaatgaacgg
cgaagtgtgtctaccactctgaagaacccgcgggagcttactgtcgaacgcactcatgtggatcgtgtgcgagcgccccctcgtgtgctccctgagatgtagctgtcttctcgtctgtgtgcccgctgctgtgagatagtagtctgtgtggatcagctgtgtgaacgcgacatt
gtgtcttctctgcgcgcatgtcagaaaaaacgaacccctcactga

More info: Please contact us for any question: contact@proteogenix.fr