

CLDN16-claudin 16 Gene

Reference: PTXBC069662

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOMG3; PCLN1; claudin-16; PCLN-1; hypomagnesemia 3, with hypercalciuria and nephrocalcinosis; paracellin-1; claudin 16		

For research use only.

GENE DESCRIPTION	claudin 16
GENE ID	10686
ORIGIN SPECIES	Human
ACCESSIONS	BC069662
VECTOR	pDONR223
NCBI SYMBOL	CLDN16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacctccaggacccactgttggttacagcctgtttgtatttcttactgcaactcaagacacctgcagcaggcgtgagaaaaagtaaaagaccagttatcttcacattgccaggtaccagaaacacagaagactgacacccgccactaaagtggggccagggtgtgtgtcccatgtt
gccatcctgatggcgtgttgcccaatgagggatcttctcaatcacgcttgccttttctctctgctgggtttttgattgtggccacctggactgactgttgatggtgaatgctgatgactctctggagggtgagcacaataatgccgagccctctgttggaatgctgcacaatgctt
ttgatgggattcgacctgtgatgagtagcattccatacttgcggagcatcccttgaagctggtgtaactcagcgtgttgatgattactgcagataattctagctgggtttggatttccacctgctccttgggtcttgactgcgtgaaattccctcctgatgagccgtacataaagtcgcatctgtt
gttgcctggagccaggttactaatagcaggtacccaggaatcattgctctgtgtgtatgctgtgtatgtgtggaactgttctactttgtttgcacaataattcttggatccaataaattggttggtcctgttggctcggatggctgggtctctgggtgttcttttgctggagctgttct
cacctgctgttatatcttttaaatggttgacctgagagaactatccttattccttgaggaaagcctattcagcgcgggtgttccatggccaagtacatactcagccctcgacagagacggccaaatgtatgctgtagacacaagggtgtaa