



Reference: PTXBC070098

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized family 31 glucosidase KIAA1161; NET37		
For research use only.			

GENE DESCRIPTION	KIAA1161
GENE ID	57462
ORIGIN SPECIES	Human
ACCESSIONS	BC070098
VECTOR	pDONR223
NCBI SYMBOL	KIAA1161
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr