

PRAP1-proline-rich acidic protein 1 Gene

Reference: PTXBC071872

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO1195; UPA; proline-rich acidic protein 1; epididymis tissue protein Li 178; proline rich acidic protein 1		

For research use only.

GENE DESCRIPTION	proline-rich acidic protein 1
GENE ID	118471
ORIGIN SPECIES	Human
ACCESSIONS	BC071872
VECTOR	pDONR223
NCBI SYMBOL	PRAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggagctcctctgtcaccagcctggtggttctgctgtgggagcaggtgcagtcaccagcaccgaagtcctatcaagatcaagtcacaactggcctcagagcaggaccagagaaggcctggggcggcctgtggtggagcctccggagaaggacgaccagct
ggtggtgctgttccctgtccagaagccgaaactctgaccaccgaggagaagccacgagtcagggcaggggccccatcctccagcaccgaagcctggatggagaccgagacacctgggcccgtgctctgagtcgccgcccaccatgacagcctgtaccacctccgctg
aggaggaccagggcaggagaggccccggtggtggtgatgccaaatcaccaggtgctcctgggaccggaggaagaccaagaccatctaccacccccagtag