



Reference: PTXBC073163

|                |                                                                                                                                                                                                                          |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                      | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | guanylate-binding protein 2; GTP-binding protein 2; guanine nucleotide-binding protein 2; guanylate binding protein 2, interferon-inducible; interferon-induced guanylate-binding protein 2; guanylate binding protein 2 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | guanylate binding protein 2, interferon-inducible                                                |
| <b>GENE ID</b>          | 2634                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC073163                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | GBP2                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]