

## C1orf212-chromosome 1 open reading frame 212 Gene

Reference: PTXBC073851

|                |                                                                                                      |                     |                |
|----------------|------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | alternative protein C1orf212; UPF0767 protein C1orf212; C1orf212; small integral membrane protein 12 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 1 open reading frame 212                                                               |
| <b>GENE ID</b>          | 113444                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC073851                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C1orf212                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggtcctactgtgggtggtgactgtctgtacctgttgacagagcccccagaccacatctaattggtttgttcttctgcccctccacaccacacagccacacaaatactggctgctcttgatggccaggcagaccagcagcagccaggggccagtgaagaggaaggc  
cgcatctgtgtgtgtggccacaagcactcaggcatctgagtttactgtgctgctgggagagagttatgagatgaacattggctgtcaatctctgtggcagcggtttggcctctagtgggaatggctgggatttggcgcttgcctttaggaggatacctgcatgtctagttccagtc  
tgcactggaaagaattcaaatatgcacctggctcccttcaatttttgcctatctttgtgctcattcttactgaaatctgtctgtcagctcaggaatgggattccccaggaaaggaaagcaccttttctgttctgggaagccagactgttcaatttgggcagggaacatgtgcctctgaat  
ttgcttga