

ING3-inhibitor of growth family, member 3 Gene

Reference: PTXBC073865

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Eaf4; ING2; MEAF4; p47ING3; inhibitor of growth protein 3; inhibitor of growth family member 3		

For research use only.

GENE DESCRIPTION	inhibitor of growth family, member 3
GENE ID	54556
ORIGIN SPECIES	Human
ACCESSIONS	BC073865
VECTOR	pDONR223
NCBI SYMBOL	ING3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgacctagaagactatctggaatgattgagcagcttcctatggatctgcgggaccgttcacggaatgcgcgagatggacctgcaggtgcagaatgcaatggatcaactagaacaagagtcagtgaaattcttatgaatgcaagaaaaataaacctgagtgagggaagagca
aatggcatccatcaaaaaagactactataaagctttggagatgcagatgagaaggttcagttggcaaacagatataatgacttgtaagtaaaagtaattgtgcatactgtgccataa