

STX16-syntaxin 16 Gene

Reference: PTXBC073876

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SYN16; syntaxin-16; syntaxin 16		
For research use only.			

GENE DESCRIPTION	syntaxin 16
GENE ID	8675
ORIGIN SPECIES	Human
ACCESSIONS	BC073876
VECTOR	pDONR223
NCBI SYMBOL	STX16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaccaggcggttaaccgacgtttctgtgttgcggaataattccatccaaaaccggcagctgttagccgagcaagtgaagtagtcacatcacctccagccctctgattcacgtagcattgctgcgttctgatgaccgtatggcactgggttcaggcatcagcttagatccagaagc
agcgtattggtgacaaaacggccacctcctaagtggtggatggagtgatgaaattcagtagatgttgccgggattaagcagaagatgaaagaattggccagccttcagacaagcatttaaacagaccaccctggatgacagcagcgaaggaacatgccattgagataactac
ccaagagatcactcagctctccacaggtgccagcgtgccgagccctgccgagccggccggcgtgctccgagcagagggggcgtgcttgggaacgtggtggcctcgtggcgcaggccctgcaggaaacttccaccagcttccggcagcagctcaggctactca
aacgcatgaagaalcgagaggaagatccagcatttttcgacacatcagtagcactaatggatgatggagacgataaacctctttaccatcgggggttttacagaggaccagttagtcttggaggagcagaacacactgatggtggaagagcgggaacgagagattccagattgtacag
tcatttctgacctgaatgaatattcagggacttagggcgatgattgtagaacagggtacagtccttgacagaattgactataacgttgaacagtcctgtatcaaaactgaagatggttgaacagcttcacagcagaacagatcaaaagaagaatcgaagatgcttgtatttaatat
tattgtcatcatcattgtgctcattgtgtcctcgttgccgtgaagtctcgtataa