



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

VPS16-vacuolar protein sorting 16 homolog (*S. cerevisiae*) Gene

Reference: PTXBC073959

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

VPS16, CORVET/HOPS core subunit; hVPS16; vacuolar protein sorting-associated protein 16 homolog; vacuolar protein sorting 16 homolog; vacuolar protein sorting protein 16

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 16 homolog (S. cerevisiae)
GENE ID	64601
ORIGIN SPECIES	Human
ACCESSIONS	BC073959
VECTOR	pDONR223
NCBI SYMBOL	VPS16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr