



Reference: PTXBC073999

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GGL; GGT 5; GGT-REL; GGTLA1; gamma-glutamyltransferase 5; gamma-glutamyl cleaving enzyme; gamma-glutamyl transpeptidase-related enzyme; gamma-glutamyl transpeptidase-related protein; gamma-glutamyltransferase-like activity 1; gamma-glutamyltranspeptidase		

For research use only.

GENE DESCRIPTION	gamma-glutamyltransferase 5
GENE ID	2687
ORIGIN SPECIES	Human
ACCESSIONS	BC073999
VECTOR	pDONR223
NCBI SYMBOL	GGT5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr