

YPEL1-yippee-like 1 (Drosophila) Gene

Reference: PTXBC074501

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FKSG3; protein yippee-like 1; DiGeorge syndrome-related protein; yippee like 1		

For research use only.

GENE DESCRIPTION	yippee-like 1 (Drosophila)
GENE ID	29799
ORIGIN SPECIES	Human
ACCESSIONS	BC074501
VECTOR	pDONR223
NCBI SYMBOL	YPEL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaaatgacaaagtcacaaacttccaagcgtatctgccgaactgtcacgaactgtacagctgtatccactgcagagcacacctggccaatcatgacgagctcatctcaagtccttcaggaggagccaggagcgcctacctcttcaattccgtggtgaacgtgggctgcggccc
tgcagaggagagggctcttctcaccggctgcatgcggtgcccacatctactgcgagaactgcaagaccacgctcgggtggaaatagagcatgcctttgagagcagtcagaaatataagggaagaaattcatcattgacctgtctcatatgatcaagacaatggctgggagtaa