

ARHGAP5-Rho GTPase activating protein 5 Gene

Reference: PTXBC075799

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GFI2; RhoGAP5; p190-B; p190BRhoGAP; rho GTPase-activating protein 5; growth factor independent 2; p100 RasGAP-associated p105 protein; p105 RhoGAP; rho-type GTPase-activating protein 5; Rho GTPase activating protein 5		

For research use only.

GENE DESCRIPTION	Rho GTPase activating protein 5
GENE ID	394
ORIGIN SPECIES	Human
ACCESSIONS	BC075799
VECTOR	pDONR223
NCBI SYMBOL	ARHGAP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtgccgccgccaccgccggggccgtgccgttgaggaggagacggaggagaccgacgttgtaggggttatgtaccgaaggactctaccgtgtcagcggggaataaaactgaccaagacaatattcaaaagcagtttgatcaagatcataatcaatctagtgcaatggaagtaa
cagtaaatgtgtagctggagcccttaaaagctttcttgcagatctgccagatcccttaattccatattctctcatccagaactattggaagcagcaaaaatccggataaaacagaacgtcttcacgttgaagaaattgtaagaatttcacctgtgaaactatgatgtattcagatacgtgata
acacatctaaacagggtagtcagcaataaaatcaactaatgacagcagacaactatccatctgttttggccaacttgatgagacctgattttgaaaatcgagagtttctgtactactaagattcatcaatctgtgtgaaacattcattcagcagtgctcagttttctttacaatggagaa
attgtagaacgacaaacattgtggctctccaccacctcaaaaccaggacagttgtggaaccaatgtgccacttcagttgccgccaccattgcaacctcagctgatacaaccacaattacaacggatcctcttggtattatga