

HCFC1R1-host cell factor C1 regulator 1 (XPO1 dependent) Gene

Reference: PTXBC080178

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HPIP; host cell factor C1 regulator 1; HCF-1 beta-propeller-interacting protein; host cell factor C1 regulator 1 (XPO1 dependent)		

For research use only.

GENE DESCRIPTION	host cell factor C1 regulator 1 (XPO1 dependent)
GENE ID	54985
ORIGIN SPECIES	Human
ACCESSIONS	BC080178
VECTOR	pDONR223
NCBI SYMBOL	HCFC1R1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcctgcagcagcccttcagcagagggcccccaggaggggcccccagcgcctcccgccggcgcccttgggggtgacttggggcctggagccaggagcctctgcgcaagcagttctgtctgaggagaacatggccaccactctctcaactcagcctgcacaatgaccaccc
tactgcagccccccatgacatttccccagccctgcccccaactcaggagcccttgccttgagctgtcttctgtgcgtatcctggcagcctcctcctgagccctcctgtctgctgaggctggggacacccccagtcctccctacccctgcaacccccagctggggacataatggagctct
ga