



Reference: PTXBC080183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein KIAA0895-like; KIAA0895 like		
For research use only.			

GENE DESCRIPTION	KIAA0895-like
GENE ID	653319
ORIGIN SPECIES	Human
ACCESSIONS	BC080183
VECTOR	pDONR223
NCBI SYMBOL	KIAA0895L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]