

OR3A1-olfactory receptor, family 3, subfamily A, member 1 Gene

Reference: PTXBC096187

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OLFRA03; OR17-40; OR17-82; OR40; olfactory receptor 3A1; olfactory receptor 17-40; olfactory receptor OR17-15; olfactory receptor family 3 subfamily A member 1		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 3, subfamily A, member 1
GENE ID	4994
ORIGIN SPECIES	Human
ACCESSIONS	BC096187
VECTOR	pDONR223
NCBI SYMBOL	OR3A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccagaatcggggccaatggaacagtcattgctgagttcatcctgctgggctgctggaggcgccagggtcagccagttgtctttgtgctctctctttgctacctgtgcaggtcaggggcaacctcagcatcctggcagctgtcttggaggcccaactccacacccc
atgtacttctctctggggaacctatcagtgctggatgttgggtgcatcagcgtcactgttccatcaatgttgagtcgtctctgctcccgcaagcgtgcagttccctgtggggcctgcttaccagctcttcttccatctgttctgttgagtgagactgcttctgctgaccgcatggcctatgac
cgattctggccatctgcccggccctcacctacagcaccgcatgagtcagacagtcagaggatgttgggtgctgctcctgggctgtgtcttccaccaacgcactgacccacactgtggccatgtccacgctcaactctgtggcccaatgtgatcaatcacttctactgtgacctccac
agctcttcagctctctgctcagcaccacccaactcaatgagctgctgttttgggttttataatggcaggtaccgccatggctctcattgtcatctctatccacgtggcagctgagctcgcgaattcgctctgtagagggcagggaagaagccttctccatgtggctccacctc
actgtggtgcatattctatgttcaggtatctttaactatgctgactgggttcaaccaagcttccagacaaggataaagctgttggaaatttcaacactgtcatcaatcccatgctgaaccaatcatctacagcttcagaacccctgatgtgcagagtgccatctggaggatgctcacaggag
gcggtcactggctga