



Reference: PTXBC096238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPYR1; NPY4-R; PP1; neuropeptide Y receptor type 4; pancreatic polypeptide receptor 1; neuropeptide Y receptor Y4		

For research use only.

GENE DESCRIPTION	pancreatic polypeptide receptor 1
GENE ID	5540
ORIGIN SPECIES	Human
ACCESSIONS	BC096238
VECTOR	pDONR223
NCBI SYMBOL	PPYR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]