



Reference: PTXBC096294

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CK-2e; K2e; KRT2A; KRT2E; KRTE; keratin, type II cytoskeletal 2 epidermal; cytokeratin-2e; epithelial keratin-2e; keratin 2, type II; keratin-2 epidermis; keratin-2e; type-II keratin Kb2; keratin 2		

For research use only.

GENE DESCRIPTION	keratin 2
GENE ID	3849
ORIGIN SPECIES	Human
ACCESSIONS	BC096294
VECTOR	pDONR223
NCBI SYMBOL	KRT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]