



Reference: PTXBC096298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACHM2; CCNC1; CCNCa; CCNCalpha; CNCG3; CNG3; cyclic nucleotide-gated cation channel alpha-3; CNG channel alpha-3; CNG-3; cone photoreceptor cGMP-gated channel alpha subunit; cone photoreceptor cGMP-gated channel subunit alpha; cyclic nucleotide gated chan		

GENE DESCRIPTION	cyclic nucleotide gated channel alpha 3
GENE ID	1261
ORIGIN SPECIES	Human
ACCESSIONS	BC096298
VECTOR	pDONR223
NCBI SYMBOL	CNGA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]