



Reference: PTXBC096312

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTN; BTN1; butyrophilin subfamily 1 member A1; bK14H9.2 (butyrophilin, subfamily 1, member A1)		

For research use only.

GENE DESCRIPTION	butyrophilin, subfamily 1, member A1
GENE ID	696
ORIGIN SPECIES	Human
ACCESSIONS	BC096312
VECTOR	pDONR223
NCBI SYMBOL	BTN1A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]