



Reference: PTXBC096324

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FEX; GPR49; GPR67; GRP49; HG38; leucine-rich repeat-containing G-protein coupled receptor 5; G-protein coupled receptor 49; G-protein coupled receptor 67; G-protein coupled receptor HG38; orphan G protein-coupled receptor HG38; leucine rich repeat contain		

For research use only.

GENE DESCRIPTION	leucine-rich repeat-containing G protein-coupled receptor 5
GENE ID	8549
ORIGIN SPECIES	Human
ACCESSIONS	BC096324
VECTOR	pDONR223
NCBI SYMBOL	LGR5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]