

HIST1H2BB-histone cluster 1, H2bb Gene

Reference: PTXBC096728

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2B.1; H2B/f; H2BFF; histone H2B type 1-B; H2B histone family, member F; histone 1, H2bb; histone H2B.1; histone H2B.f; histone cluster 1, H2bb; histone cluster 1 H2B family member b		

For research use only.

GENE DESCRIPTION	histone cluster 1, H2bb
GENE ID	3018
ORIGIN SPECIES	Human
ACCESSIONS	BC096728
VECTOR	pDONR223
NCBI SYMBOL	HIST1H2BB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgaacctctaagtctgtccagccctaaaaagggttctaagaaggctactaaggcgcagagaaggatgtaagaagcgttaagcgcagccgcaaggagagctattctatctatgtgtacaagggtctgaagcaggtccacccgacaccggcatctccaaggccatg
gggatcatgaattcctctgtaacgacatcttcgagcgcacgcggggcagggcttctgcctggctactacaataagcgtcgaccatcacctccaggagattcagacggctgtgcgcctgctgctgcctggggagctggctaagcatgctgtgtccgagggcactaaggcagttacca
agtacactagctctaaataa