



Reference: PTXBC098347

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCDH-BETA2; protocadherin beta-2; protocadherin beta 2		
For research use only.			

GENE DESCRIPTION	protocadherin beta 2
GENE ID	56133
ORIGIN SPECIES	Human
ACCESSIONS	BC098347
VECTOR	pDONR223
NCBI SYMBOL	PCDHB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr