



Reference: PTXBC099634

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AKR6A9; KCNA3.1B; KV-BETA-3; voltage-gated potassium channel subunit beta-3; K(+) channel subunit beta-3; potassium channel, voltage gated subfamily A regulatory beta subunit 3; potassium channel, voltage-dependent, beta-3 subunit; potassium voltage-gated		

For research use only.

GENE DESCRIPTION	potassium voltage-gated channel, shaker-related subfamily, beta member 3
GENE ID	9196
ORIGIN SPECIES	Human
ACCESSIONS	BC099634
VECTOR	pDONR223
NCBI SYMBOL	KCNAB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]