



Reference: PTXBC099652

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GluD2; SCAR18; glutamate receptor ionotropic, delta-2; gluR delta-2 subunit; glutamate receptor delta-2 subunit; glutamate receptor, ionotropic, delta 2; glutamate ionotropic receptor delta type subunit 2		

For research use only.

GENE DESCRIPTION	glutamate receptor, ionotropic, delta 2
GENE ID	2895
ORIGIN SPECIES	Human
ACCESSIONS	BC099652
VECTOR	pDONR223
NCBI SYMBOL	GRID2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]