

HIST3H2BB-histone cluster 3, H2bb Gene

Reference: PTXBC100855

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2Bb; histone H2B type 3-B; H2B type 12; histone 3, H2bb; histone cluster 3, H2bb; histone cluster 3 H2B family member b		

For research use only.

GENE DESCRIPTION	histone cluster 3, H2bb
GENE ID	128312
ORIGIN SPECIES	Human
ACCESSIONS	BC100855
VECTOR	pDONR223
NCBI SYMBOL	HIST3H2BB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagaccgtccaaatcgctcctgcgcaccaagaagggttctaaaaaggctgtcaccaggcacagaagaaggacggcaagaagcgccgcaaggagagctattctatctacgtgtacaagggtgtgaagcaggtgcaccccgacaccggcatctcgtccaaggc
catgggcatcatgaactcctctgtaaatgacatctcgagcgcatcgccagcgaggcctcccgctggcacactacaacaagcgctccaccatcacgtcccgcgaagtgcagacggcgttcgctgctgctgcccggcgagctgccaagcacgcccgttccgagggcaccaaggct
gtcaccaagtacaccagctccaagtga