

MGC119295-similar to Williams-Beuren syndrome critical region protein 19 Gene

Reference: PTXBC100973

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to Williams-Beuren syndrome critical region protein 19
GENE ID	441273
ORIGIN SPECIES	Human
ACCESSIONS	BC100973
VECTOR	pDONR223
NCBI SYMBOL	MGC119295
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtggggacgaatctgaggagtcgtggaggagccacggaaggtctcgccttgagcctgaggagatctgggtggcggagatctgtgtgacctcaagatgaagctgaagcgacggcgagtgctcgtcgtcctgagcaccagggcctcaacaggtgctt
gaggatcctgtcattaaaagattcttggcctgggacaagatctgagggtgtcggacaagtatctcctggtatggtcatagcgtatttcagccggcggtctccctcctggcaataccaacgcattcatttctcctggctctctacctggccaatgacatggaggagcagcaggact
ccaaacaaacatcttccactctgtataggagaacgctctcgcataccttgcctcgttaagcgttgggtccagttaggccattccatgaacccgaggccaggagaacccgctctcgcataccttgcctcgttaagcgtcgggtccagttataccgttccacgaacccgaggccagga
agaaccgctctcgcataccttgcctcgttaagcgtcgggtccagttatacgttccatgaactcaggggccaggagaacccgctctcagatagtcctgttccagaacgacggttccacttctctgttccatgagctgcagggttgggttccccagaggagttggagagatccaggtta
tgaccagagacgactgggtgtggcgcgagatcgcgtcaccttccatg