

HTR3D-5-hydroxytryptamine (serotonin) receptor 3 family member D Gene

Reference: PTXBC101090

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	5HT3D; 5-hydroxytryptamine receptor 3D; 5-hydroxytryptamine (serotonin) receptor 3 family member D; 5-hydroxytryptamine (serotonin) receptor 3D, ionotropic; 5-hydroxytryptamine receptor 3 subunit D; serotonin 5-HT-3D receptor; serotonin receptor 3D		

For research use only.

GENE DESCRIPTION	5-hydroxytryptamine (serotonin) receptor 3 family member D
GENE ID	200909
ORIGIN SPECIES	Human
ACCESSIONS	BC101090
VECTOR	pDONR223
NCBI SYMBOL	HTR3D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggccatcaggcgaggtgcaggccagccctacgtggtaactttctgtgccagtggtcattctgattgccatcgatgccctcagttctacctgccactggaaagtgggaattgtgccccattcaagatgactgttctgctgggtacagcgtctctctgctcatgatgaatgacttg
ctccagccactagcaactcatcatgtttcactagtagctctcatcatcaagagacaaaagcgagggtgtctacttgcctgtgctctgctgatgtgtggcgagcctgtggagaccattctcatcaccacgtgctgcagtggtgcccaccacccagccctacctctgctcggtgg
ctccactcctgctgctgcactgcaccggccaaggagatgctgtccactgcgccccagagggaataaaggccgggtctcaccacccacacgtgcccgtgtgaaggagccagaggtatcagcagggcagatgccaggccctggggaggcagagctgacagggggctca
gaatggacaaggccagcgggaacacgaggccagagcagcactcgggtgagctgtgggtgcagttcagccacgcgatggagccctgctcttccgctctacgtgctctcatggcctctccatcatcaccgtcatatgcctctggaacacctag