



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GPR120-G protein-coupled receptor 120 Gene

Reference: PTXBC101175

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

GPR120; BMIQ10; GPR129; GT01; O3FAR1; PGR4; free fatty acid receptor 4; G-protein coupled receptor 120; G-protein coupled receptor 129; G-protein coupled receptor GT01; G-protein coupled receptor PGR4; omega-3 fatty acid receptor 1

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 120
GENE ID	338557
ORIGIN SPECIES	Human
ACCESSIONS	BC101175
VECTOR	pDONR223
NCBI SYMBOL	GPR120
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr